LUXAMED®

LUXAMED® Untersuchungsleuchte
LUXAMED® Examination Lamp
Éclairage d'examen LUXAMED® LED

| | |
|---------------------|----|
| Deutsch | 2 |
| English..... | 6 |
| Français..... | 10 |
| EMC Guidelines..... | 14 |

LUXAMED GmbH & Co. KG * Daniel-Weil-Str.3 * 89143 Blaubeuren * Germany
Tel.: +49(0)7344 92905 0 * Fax: +49(0) 7344 92905 10
E-Mail: info@luxamed.de * www.luxamed.de

- DEUTSCH -

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Minimaler Biegeradius: 60 mm

Kaufdatum

Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

| | Separation distance according to frequency of transmitter in m | | |
|------|--|--|---|
| | 150 kHz to 80 MHz $d = \{ 3,5/\sqrt{1} \} \sqrt{P}$ | 80 MHz to 800 MHz $d = \{ 3,5/E1 \} \sqrt{P}$ | 800 MHz to 2,5 GHz $d = \{ 7/E1 \} \sqrt{P}$ |
| 0,01 | 0,1 | 0,1 | 0,2 |
| 0,1 | 0,4 | 0,4 | 0,7 |
| 1 | 1,2 | 1,2 | 2,3 |
| 10 | 3,7 | 3,7 | 7,4 |
| 100 | 11,7 | 11,7 | 23,3 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1—At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2—These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

| Guidance and manufacturer's declaration-electromagnetic immunity | | | |
|---|---|---|--|
| The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an electromagnetic environment. | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment-guidance |
| Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3 | 3 Veff 150 kHz to 80 MHz 3 V/m 80 MHz to 2,5 GHz | V1 = 10 Vemk @ 150 KHz to 80 Mhz E1 = 10 V/m @ 80 MHz to 2,5 GHz | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \{ 3,5/\sqrt{V1} \} \sqrt{P}$ $d = \{ 3,5/E1 \} \sqrt{P}$ 80MHz to 800 MHz $d = \{ 7/E1 \} \sqrt{P}$ 800 MHz to 2,5 GHz 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, c should be less than the compliance level in each frequency ranged Interference may occur in the vicinity of equipment marked with the following symbol:  |
| NOTE 1 - At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people. a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m. | | | |

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- ENGLISH -

 To ensure proper use of this LUXAMED® product, please read this user manual carefully and retain for future reference.

LUXAMED® Examination Lamp

Purpose:

The LUXAMED® examination lamp is only to be used in medical practice. This examination lamp is operated by local electricity grid and is used for general auxiliary illuminating the examination field. The examination lamp is only to be used by qualified medical personnel. Use the device in accordance with its designated use and for its intended purpose only.


Warning:

1. Possibly hazardous optical radiation emitted from this product. Do not stare at operating lamp. May be harmful to the eyes. Do not direct the examination lamp into eyes.
2. Please note that the head of the examination lamp may heat up after longer use.
3. This examination lamp may only be used, if dangers can be excluded in the event of a possible failure. Another main illumination must be presented.
4. This examination lamp may not be used in MRT, areas with increased magnetic field and potentially explosive areas.
5. Do not cover the head of the examination lamp.
6. Near flammable anaesthetics (class AP) or flammable anaesthetics with oxidizing agents (class APG), the examination lamp may not be used.
7. Do not bend or tension power cable. Further, please make sure that the power cable does not present a tripping hazard.
8. Please make sure that the examination lamp is not contaminated during operation.
9. Make sure that there is no tripping hazard caused by the electric line.
10. When switching on the lamp, a short, single flash may occur, due to integrated current monitoring.
11. To avoid endangering a patient, do not perform testing or maintenance when the machine is in use.
12. MAINS plug is used as the disconnect device, which shall remain readily operable.
13. The product is indoor use only, can only be in used in overvoltage II and less than 3000m.
14. Do not bend the metal hose too much, never less than a bending radius of 60 mm.
15. All serious incidents that have occurred with the product must be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established

Commissioning:

Package contents:

1 x base

1 x stand tube

5 x base tube

| Guidance and manufacturer's declaration-electromagnetic immunity | | | |
|---|---|---------------------------------|---|
| The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment. | | | |
| Immunity test | IEC 60601 - test level | Compliance level | Electromagnetic environment-guidance |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±6 kV contact ±8 kV air | ±6 kV contact ±8 kV air | Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Electrical fast transient/burst IEC 61000-4-4 | ± 2 kV for power supplies ± 1 kV for in- and output lines | ± 1 kV for in- and output lines | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ± 1 kV differential mode ± 2 kV common mode | n/a | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 s | n/a | A typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 3 A/m | 3 A/m | The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low. |
| NOTE - UT is the a.c. mains voltage prior to application of the test level. | | | |

| Guidance and manufacturer's declaration-electromagnetic emissions | | |
|---|------------|---|
| The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment. | | |
| Emissions test | Compliance | Electromagnetic environment-guidance |
| RF emissions CISPR 11 | Group 1 | The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Complies | |

Assembly of the stand (from left to right):

First, mount the stand by connecting the stand tube with plastic foot part. Then, put the metal plate with the bulges in direction to stand tube and insert screw. Tighten the screw just enough to insert base tubes sideways. In this process, the elevations have to be protruded into the hole of the base tubes.

After fitting all 5 base tubes, tighten the screws firmly until the stand is placed in a stable position. Press the arm of the examination lamp with the snap-in into the stand. The lamp must be fixed to the stand by the lateral stud screw. To do this, turn the stud screw clockwise with the supplied Allen key.

The examination lamp can now be connected with the local power supply. By pressing the main power switch, the examination lamp can either be switched on or switched off. Please ensure to maintain a minimum distance of 200 mm to the examination field.

Cleaning:

Before cleaning, turn the examination lamp off and pull the plug.
Clean with a soft, lint-free cloth, slightly moistened with alcohol, if required. Do not touch electrical contacts.
Always clean the examination lamp after use.

Disinfection:

Gemäß der Empfehlung der Kommission für Krankenhaushygiene und Infektionsprävention (KRINKO) des According to recommendations by the Commission for Hospital Hygiene and Infection Prevention (KRINKO) at the Robert Koch Institute (RKI) and the Institute for Drugs and Medical Devices (BfArM):
Only use disinfectants with proven antibacterial, antiviral and fungicidal effectiveness.
Bundesgesundheitsblatt 2012 · 55:1244-1310; DOI 10.1007/s00103-012-1548-6.
In principle, CE-marked alcohol based instrument disinfectants (2-propanol) are suitable.
Please read the operating instructions of the used disinfectant carefully.

Wipe the head of the examination lamp with a disinfectant impregnated cloth. Comply with the exposure time according to the operating instructions of the disinfectant. Then wipe possible residues.

Sterilisation:

Not approved for sterilisation.

Maintenance:

No maintenance is required. In case of damage, the LED examination lamp can be repaired by the manufacturer.

Inspection and functional testing:

Check the examination lamp and the power cable carefully before use. Switch on the examination lamp in order to guarantee its proper function. Immediately after switch on the lamp there may be a single flickering of the light which does not impair the function and safety. This short and single flickering is due to current spikes at turning-

on which get filtered by the AC adapter. This effect protects the LEDs and increases the life time of the LED light source..

Disposal:

The arm of the examination lamp must be disposed of at a separate collection of electric and electronic equipment.
 The stand can be disposed of with metal scrap.

Storage, transport and operating conditions:

Ambient temperature between -15° and +45° C.
 Rel. humidity between 20% and 95%.

Technical data:

| | |
|-------------------------|---------------------------------|
| REF U1.111.014 | Examination Lamp LED with stand |
| REF U1.111.114 | Examination Lamp LED |
| Weight: | 6,2 kg (with stand) |
| Mains connection: | 100-240 V AC, 50/60 Hz, 0.6 A |
| Lamp life: | approx. 50.000 hours |
| Product life: | ≥ 5 years |
| Illuminance: | 50.000 Lux |
| Colour temperature: | approx. 4.000 K |
| Operating temperature: | + 5°C until + 40°C |
| Rel. humidity: | 10% until 75% |
| Air pressure: | 700 hPa until 1070 hPa |
| Working height max.: | 1650 mm |
| IP protecting class: | IP 20 |
| Length of flexible arm: | 900 mm |

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60 mm min.

Symbol key

| | |
|--|--|
| | Manufacturer |
| | Date of manufacture |
| | Attention to the user manual |
| | Caution |
| | Separate collection of electric and electronic scrap (European WEEE Directive) |
| | In Conformity with EU guidelines |
| | This is a medical device |
| | Interference may occur in the vicinity of equipment marked with this symbol. |

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